



Be Right™



Nitrate test kit, model NI-11

Продукт №: 146803-EU
EUR Цена: Свържете се с нас
Не е показана дата на доставка

Опасен



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Спецификации

Parameter:	Nitrogen, Nitrate
Smallest increment:	1.0
Бр. тестове:	100
Вид на корпуса:	B
Измерване, диапазон:	Single Range
Метод, име:	Color disc/ Cadmium reduction
Модел:	NI-11
Модел на тестов комплект:	NI-11**
Параметър:	Nitrogen, Nitrate - As NO_3^- -N
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 40 mg/L NO_3^- -N